

Easy Rider Compact Small

Compact tricycle for children

NEW



Van Raam's compact tricycle

The Easy Rider Compact Small is a compact tricycle with a low step-through frame. The low center of gravity provides great stability while cycling. The seat with backrest offers extra riding comfort. The bike allows for a forward pedaling motion, enabling you to apply more power when pedaling.

The Easy Rider Compact Small can be easily converted into an Easy Rider Compact. This allows the bike to grow with the rider.

For whom?

The Easy Rider Compact Small is suitable for children who have difficulty with stability. Individuals with an inner leg length of approximately 52 cm to 72 cm can fit on this tricycle.

”

A major advantage is that the Easy Rider Compact Small is a grow-along bike and can be converted into the Easy Rider Compact.





Benefits Easy Rider Compact Small

- ✓ Easily convertible to an Easy Rider Compact model
- ✓ Very compact tricycle with a small turning radius
- ✓ Low step-through frame for easy mounting
- ✓ Rides easier than a traditional tricycle
- ✓ Less strain on the back, neck, shoulders, and wrists

Characteristics Easy Rider Compact Small

- With pannier rack suitable for single-sided panniers
- Comfortable seat with adjustable backrest
- Sporty and modern appearance
- A stable, low seating position ensures a safe feeling
- Fits through a standard door opening

Options and pedal support

With various options, the bike can be customized to your preferences. For example, think of foot fixation, a fixation vest, or pedal adjustments.

Do you want an electric tricycle? The Easy Rider Compact Small can also be equipped with pedal assistance. The maximum supported speed for this bike is 15 km per hour.

Technical specifications

Entry height (cm):	32
Bike weight ¹ (kg):	37,5
Max. weight pannier rack (kg):	10
Total bike length (cm):	173
Total bike width (cm):	73
Average range ² (km):	61
Max. user weight (kg):	80
Inner leg length (cm):	± 52-72

1. Without the electrical system

2. With default battery



More information?

Do you want more information about the Easy Rider Compact Small tricycle? Then please contact Van Raam.

☎ +31 (0)315 - 25 73 70

✉ info@vanraam.com

🌐 www.vanraam.com